

High-Throughput Cloning Vectors (including MBP/T4L Mutants) Multiple Cloning Sites

**Note: All high-throughput cloning vectors listed in this document are AMPICILLIN resistant

*All cloning in this vectors can be done with Froward BamHI and Reverse XhoI or STOP-XhoI

HT1: N-His-TEV-Protein-STOP Original Base Vector was pet16b

*Sequence with T7 & T7Term Primers

Bgl II T7 Promoter lac operator XbaI
AGATCTCGATCCCGCAAATTAATACGACTCACTATAGGGGAATTGTGAGCGGATAACAATTCCCCTCTAGAAATAATTTTCTTTAAGTT

RBS NcoI His Tag TEV Cleavage Site ↓ BamHI
TAAGAAGGAGATATACCATGGGCCGGGGCAGCCATCATCATCATCATGGCGCGCAAACCTGTATTTTCAGGGCCCGGGATCC
MetGlyArgGlySerHisHisHisHisHisHisHisGlyGlyGluAsnLeuTyrPheGlnGlyProGlySer

EcoRI XhoI *Note the XhoI site is out of frame – you need to add a STOP codon
CCGAATTCCCGGGTTCGACTCGAGCGGCCGC...
ProAsnSerArgValAspSerSerGlyArg

HT2: Protein-His-STOP Original Base Vector was pet21b – T7 tag removed from N-terminus

*Sequence with T7 & T7Term Primers

Bgl II T7 Promoter lac operator XbaI
AGACTCCGATCCCGCGAAATTAATACGACTCACTATAGGGGAATTGTGAGCGGATAACAATTCCCCTCTAGAAATAATTTTGTTTAACT

RBS NdeI BamHI EcoRI XhoI His Tag
TTAAGAAGGAGATATACATATG GGATCC GAATTC CTCGAG CACCACCACCACCACCCTGA
Met GlySer GluPhe LeuGlu HisHisHisHisHisHisHis *

HT4: N-Trx-His-TEV-Protein – TEV site replacing Enterokinase site in pet32b

*Sequence with S•Tag primer #69945-3 & T7 Term

T7 Promoter→ lac operator XbaI RBS
TTATACGACTCACTATAGGGGAATTGTGAGCGGATAACAATTCCCCTCTAGAAATAATTTTGTTTAACTTTAAGAAGGAGATATAACAT

Trx Tag MscI His Tag
ATGAGC ...315bp/105aa... CTGGCCGGTTCTGGTTCTGGCCATATGCACCATCATCATCATCATTCTTCTGGT
MetSer LeuAlaGlySerGlySerGlyHisMetHisHisHisHisHisHisHisSerSerGly

Thrombin ↓ S Tag Bgl II
CTGGTGCCACGCGGTTCTGGTATGAAAGAAACCGCTGCTGCTAAATTCGAACGCCAGCACATGGACAGCCCAGATCTG
LeuValProArgGlySerGlyMetLysGluThrAlaAlaAlaLysPheGluArgGlnHisMetAspSerProAspLeu

TEV Cleavage ↓ BamHI EcoRI SacI Sal I HindIII NotI XhoI
GAAAACCTGTATTTTCAGGGCCCGGGATCC GAATTC GAGCTCCGTCGAC AAGCTT GCGGCCGCACTCGAG
GluAsnLeuTyrPheGlnGlyProGlySer GluPhe GluLeuArgArg GlnAla CysGlyArgThrArg...

His Tag (Not in frame) *Note XhoI not in frame – Add STOP to reverse primer to avoid extra residues
CACCACCACCACCACCCTGA

