

## Ligation (sticky-ends)

### Special Reagents:

Crystal violet (10 mg/ml solution in water)

Agarose gel reagents (agarose, TBE)

T4 ligase (and included buffers)

Restriction enzymes and buffers

1. Prepare a standard miniprep of plasmid DNA (including target vector and plasmid containing fragment to be subcloned in to target vector).
2. **Digest 5  $\mu$ L of each miniprep in a 20  $\mu$ L reaction.** This usually means 1  $\mu$ L of each of the appropriate restriction enzyme(s), 2  $\mu$ L of 10X buffer, and 11  $\mu$ L of ddH<sub>2</sub>O.
3. Run products of this reaction out on an **agarose gel stained with crystal violet**. To prepare gel, mix and boil gel as usual, but before pouring add 10 ug/ml crystal violet (ie, 50  $\mu$ L of 10 mg/ml crystal violet solution in 50 mL gel solution). **Gel must be run in TBE with crystal violet as well!!!**
4. Visualize gel on a standard light box (no UV) and (with a razor blade) cut out bands corresponding to target vector and fragment, and transfer each to microfuge tubes.
5. Gel-purify DNA from these bands using the Qiaquick gel extraction kit (follow kit instructions exactly as written).
6. Set up ligations as follows (use 0.2  $\mu$ L PCR tubes):
  - 5  $\mu$ L vector DNA (from step 5)
  - 9  $\mu$ L insert fragment DNA (from step 5)
  - 3  $\mu$ L ddH<sub>2</sub>O (sterile)
  - 2  $\mu$ L 10X T4 ligase buffer (well-dissolved—this stuff precipitates during storage and must be vigorously mixed to re-dissolve the ATP, DTT, etc.)
  - 1  $\mu$ L T4 ligase (roughly 2.5 Weiss units)
7. Incubate at 16<sup>0</sup>C overnight (in PCR machine), or 2 hours at room temperature.
8. Heat-inactivate digestion mix (65<sup>0</sup>C for 15 minutes).
9. Transform entire digestion mix in to MC1061 heat-shock competent E. coli. Plate on appropriate selection media.